


Item information for the instrument set

**Zestaw Promocyjny
NEVADENT SP. Z O.O.**



Description:

8 rotating BUSCH instruments, disinfected, sterilisable and sorted on
1 printed bur stand BUSCH-STERI-SAFEwave 5773-01 made of medical plastic disinfected, sterilisable with red silicone band



#	Reference number	Description	Maximum rotation speed:	Labelling
023556	1SXM 021 CARBIDE RA	Carbide instrument with x-cut and slender neck	<3500 RPM	CE 0124 Co, CAS:7440-48-4
023606	1SXM 016 CARBID RA-L	Carbide instrument with x-cut, slender neck, and long shank	<4500 RPM	CE 0124 Co, CAS:7440-48-4
032130	CD 39 016 COOL-DIA FG	COOL DIAMOND with asymmetric intersecting grooves for transport more cooling medium	<450000 RPM	CE 0124
033160	MI 889 007 DIAMOND FG	Diamond instrument with small working part diameter and long, thin neck for minimal invasive preparation	<160000 RPM	CE 0124
045584	9591 140 COMP POL DIA M	First step of 2-step diamond interspersed polishing system for composite (also high filled)	<20000 RPM	CE 0124
045524	9581 140 COMP POL DIA F	Second step of 2-step diamond interspersed polishing system for composite (also high filled)	<20000 RPM	CE 0124
045526	9583 055 COMP POL DIA M	Second step of 2-step diamond interspersed polishing system for composite (also high filled)	<20000 RPM	CE 0124
045594	9593 055 COMP POL DIA F	First step of 2-step diamond interspersed polishing system for composite (also high filled)	<20000 RPM	CE 0124

General:

Safety and hygiene recommendations for rotary instruments: see back page.
At www.busch.eu this information is always displayed up to date.

Contact:

BUSCH & CO. GmbH & Co.KG
Unterkaltenbach 17-27 · 51766 Engelskirchen - Germany
Phone +49 2263 860 · Fax +49 2263 207741
www.busch.eu / mail@busch.eu